



Microvial with bottom fused integrated insert (BaseBond)

Microsampling Vials with bottom fused insert CG-Bond

Our “**CG-Bond**” micro autosampler vials are state-of-the-art sampling vials with a stable bottom shape that can be placed safely in the autosampler as well as on the laboratory bench without the use of adapters.

A classic 32x11.6mm autosampler design, the CG-Bond vials can be used universally on all common autosamplers.

They are always used when reduced sample volumes are available and a high recovery with low residual volume is needed.



Order information:

item	description	vol. (ml)
KG 09 0682	Short Thread Vial ND9, clear, 1. hydrol. class	0,3
KG 09 0683	Short Thread Vial ND9, amber, 1. hydrol. class	0,3
RR 11 0680	Crimp Neck Vial ND11, clear, 1. hydrol. class	0,3
RR 11 0681	Crimp Neck Vial ND11, amber, 1. hydrol. class	0,3
SR 11 0683	Snap Ring Vial ND11, clear, 1. hydrol. class	0,3
SR 11 0684	Snap Ring Vial ND11, amber, 1. hydrol. class	0,3

CHROMA GLOBE GbR

Schneidhausen 2

52372 Kreuzau

Tel.: 0 24 22/ 90 41 38-0

Fax: 0 24 22/ 90 41 38-7

e-mail: service@chroma-globe.de

Homepage: www.chroma-globe.de